



DIGITAL COLONY COUNTER MODEL NO. ALE - 552/553

Digital & Microprocessor Colony Counter is designed for quick and accurate counting of bacterial and mould colonies in petri dishes.

It is designed for rapid and accurate counting of bacterial and mould colonies. Simply place the petri dish on the illuminated pad and touch the dish with the pen provided to mark each colony in turn. This causes a count to be registered on the digital display and an audible tone confirms each count made. Marking the dish with the pen avoids missing colonies or double counting. The digital count on the display can be reset manually any time by pressing reset key provided.

Optimum viewing of colonies is aided by peripheral glare free illuminations. An Integral magnifying glass provides for easier counting of small colonies.

Features

- Rugged and reliable
- Easy to operate & User friendly
- Readout 0 - 9999 (Model No-552)
- Readout 0 - 999999 (Model No-553)
- Auto marker pen
- Uniform glare free illumination
- Audible confirmation of each count
- Wolfhugel glass grid with focusing facility
- 110 mm dia magnifying glass.

Accessories

- Marking Pen
- Magnifying Lens
- Adaptor
- Operation Manual
- Dust Cover

SPECIFICATION

★ Five Stages of Quality Checking Before Dispatch.

Model	ALE-552	ALE-553
Digital Display	4 Digits	6 Digits
Max. Counts	9999 max. Count.	999999 max. Count.
Dish Size	110 mm	110 mm
Magnification	1.7x	1.7x
Power	30 V, $\pm 10\%$ AC, 50 Hz, 40W	30 V, $\pm 10\%$ AC, 50 Hz, 40W
Dimensions	260 x 220 x 168 mm (L x B x H) (Approx.)	260 x 220 x 168 mm (L x B x H) (Approx.)
Weight	3 Kg (Approx.)	3 Kg (Approx.)

Note: Technical Specifications, Appearance & Model number are subject to change without prior notice.